


SCANX COMPACT


Medical & Beauty Skin Test & Analyser


- » 5" LCD touch screen
- » Integrated observation, photographing, report generation, and direct-to-printer function
- » Easy to take pictures of pathological observations and timely trace the development of the disease and treatment
- » Suitable for use in the medical, beauty and vet industry
- » Large storage capacity for images
- » Portable, lightweight, easy to carry



MM MEDICAL
AESTHETICS

 (07) 5577 9881 | 0414 970 019

 mike@beautyequipment.com.au







 www.mmmedicalaesthetics.com.au

What is a Medical & Beauty Skin Test & Analyser

Also known as a Wood's lamp, the Scanx examination is a procedure that uses transillumination (light) to detect bacterial or fungal skin infections. It also can detect skin pigment disorders such as vitiligo and other skin irregularities. This procedure can also be used to determine corneal abrasion (scratch) on the surface of your eye. This test is also known as the black light test or the ultraviolet light test.



Colour Guide

	Pale Blue: Normal & healthy skin		Brown: Pigmentation & dark spots
	White: Dead skin cells		Yellow: Oily areas of the face*
	Dark Blue: Thinner, dehydrated skin		Dark Pink or Orange: Congested pores & comedones*

* Yellow, orange or dark pink will often show as small dots (or pinpricks of light) on the face

